

Field 3D PTV v1.0

Flow Field Imaging Laboratory, University of Minnesota

Flow Field Imaging Lab



Outlines

- Product description
 - Hardware
 - Hardware overview
 - Detailed description of key components
 - Software
 - Software overview
 - Detailed description of software used
- Image Acquisition procedure
- Data processing procedure
- Demonstration cases
 - Case I: Confetti Settling
 - Case II: Snow Settling
 - Case III: Pollen Dispersion




Product Description

Flow Field Imaging Lab



Field 3D PTV System v1.0 Overview




- Hardware: a WiFi-controlled four-camera multi-view imaging system for data recording and a drone with two-color LED spaced in a fixed distance for calibration
- Software: Remote camera control for field data acquisition, EasyWand code for calibration, and OpenLPT for 3D particle reconstruction

Hardware




Hardware Overview



- Main hardware components
 - Imaging module for acquiring and storing PTV data
 - Central computer for controlling camera alignment and data acquisition
 - Calibration module for camera alignment and 3D image reconstruction
 - Lighting module for enhancing the signal in the desired sample volume (optional)
- Entire imaging hardware component can be transferred through a customized mounting platform

- Hardware components: ~\Hardware\Open3DPTV_v1p0_component list.xlsx
- Imaging module (the four boxed units): camera, motorized stage, data processing, data storage, tripods
- Lighting module: optional depending on the background; some cases (e.g., pollen dispersion) use sun light.
- Calibration module: red and green LED attached to a carbon tube, with 34.5 cm gap

Imaging Module



- Include four imaging units, protective wooden boxes for easy setup, each mounted on a tripod
- Cameras are synchronized via sync cables, with one master imaging unit controlling three slave units
- Each imaging unit includes a camera with a lens attached, a two-axis motorized stage for camera mounting and alignment, and GPS sensor. Data acquisition board are inside the wooden box enclosure




Flow Field Imaging Lab



Notes:

- Camera: 2.8 Megapixel at 95 FPS, Lens: 16 mm f/1.4
- Two-axis motorized stage: one axis for pitch control with acceleration sensor attached and one axis for yaw control with magnetometer attached
- Data acquisition board: Nvidia Jetson Nano
- Box enclosure: 20 cm x 20 cm x 20 cm in dimension (design files in ~/Hardware/2nd Gen PTV Box)
- Platform dimension: 92 x 64 cm
- design files in ~/Hardware/Gimbal system/Printed Gimbal

Imaging Module: Camera and Motorized Stage Control




- Camera, WIFI adapter, and SSD hard drive connected to Jetson Nano directly through USB ports
- Sensors (magnetometers and accelerometers) and motors (pitch and yaw) connected to a custom PCB board attached to the Jetson Nano, with header pins and wired connections designed for easy replacement
- SSD hard drive held within the boxes and stack with the Jetson Nano and board

Notes:

- PCB design file location: ~/Hardware/Gimbal system/PCB Files
- Motorized stage control includes sensors (magnetometer, accelerometer) and motors (pitch and yaw)
- Motorized stage connected to Jetson Nano through custom PCB board
- Camera (blue colored cable), SSD, WIFI through USB port

Central Computer for Data Acquisition Control

- Filed laptop for deployment
- Required Software: *Remote Camera Control* (latest version: v1.66)
 - Calculates camera spacing based on measurement volume and height
 - Integrates with the GPS control plugin to provide feedback control for motorized stage positioning
 - Manages imaging settings and camera recording functions






Notes:

- Remote Camera Control package included in ~/Software, downloadable from GitHub
- GitHub link: <https://github.umn.edu/HongFlowFieldImagingLab/remote-camera-control>

Calibration module

- Calibration Overview
 - Calibration establishes the geometric relationship between multiple cameras
 - Uses a calibration bar equipped with fixed-distance red and green lights
- User Tips
 - A drone is used to carry the bar; depending on the application, manual handling may be an option
 - Move the bar throughout the full sampling volume to ensure accurate calibration
 - Keep the bar horizontal and ensure the lights face downward during calibration



Flow Field Imaging Lab



Notes:

- Calibration bar: the calibration scale depends on the drone used. The calibration bar used in our experiment is 34.5 cm long.

Software


Flow Field Imaging Lab



Notes:

- Location: ~/Software

Software Overview



Notes:


- Step I (3DPTV system setup parameter calculator): ~/Software/3DPTV_Image_Acquisition/system setup parameter calculator
- Step II (3DPTV camera setup control): ~/Software/3DPTV_Image_Acquisition/camera setup control
- Step III (calibration image preprocessing): ~/Software/3DPTV_Processing_Pipeline/calibration image pre-processing
- Step IV (3D calibration matrix generation using easywand): ~/Software/3DPTV_Processing_Pipeline/easyWand5
- Step V (3D PTV tracking using OpenLPT): ~/Software/3DPTV_Processing_Pipeline/OpenLPT - ShakeTheBox
- Step VI (3D trajectory visualization): ~/Software/3DPTV_Processing_Pipeline/3D trajectory visualization

Step I: 3D PTV system setup parameter calculator


- MATLAB code for determining camera location and orientation

- User inputs
 - Sampling volume dimensions and location (center height to the ground)
 - Radius of cameras
 - Camera parameters
- Spherical working distance derived from sampling volume dimension
- Camera position and orientation (yaw and pitch) calculated based on above inputs


Planned camera locations



Planned camera yaw and pitch





Spherical working distance approximation



- Input parameter examples: 6 m x 6 m x 6 m with center of sample volume at 12 m above the ground.
- Camera parameters can be found in the hardware section
- Software location: ~\Software\Multicamera System Calculator

Step II: 3D PTV camera setup control

- Remote camera control GUI via SSH
- Panel for camera settings
 - Pitch and yaw alignment adjusted with feedback control GUI
 - Imaging settings: frame rate, exposure, etc.
 - Multi-camera preview streaming
- Panel for video recording
 - Simple trigger with camera synchronization
 - Recording progress display



Flow Field Imaging Lab




Notes:

- Launch the full remote camera control software and select from plugins
- Record using .raw file for high fps
- RGB mode for calibration imaging, BW for data acquisition
- Software location: ~\Software\camera setup control


Step III: Calibration image pre-processing

- Python script for pre-processing calibration images
 - Load raw images
 - Extract marker coordinates from calibration images
 - Automatically detect red and green dots
 - Save marker coordinates as .csv files
 - Export an .avi file with markers overlaid for result review

Raw calibration images



Markers identified




Notes:

- .raw format is used when recording to reach higher fps, needs decoding
- Color setting: Calibration images are RGB for identifying red and green dots, tracking images are monochrome
- Code: ~/Software/3DPTV_Processing_Pipeline/easyWand5/calibration image pre-processing

Step IV: Calibration (easyWand5)


- Wand-based multi-camera calibration tool
 - GUI for calibrating 2 or more cameras using a calibration wand
 - Based on the MATLAB Computer Vision toolbox
- Converts calibration results for use in OpenLPT
 - Recalculates calibration using the Tsai method employed in OpenLPT
 - Generates necessary configuration files for running OpenLPT



- easyWand does not support Mac with Apple Silicon
- Downloadable from: <https://biomech.web.unc.edu/wand-calibration-tools/>
- easyWand5 location: ~/Software/3DPTV_Processing_Pipeline/easyWand5
- Converting results: ~/Software/3DPTV_Processing_Pipeline/easyWand5/easyWand to OpenLPT

Step V and VI: 3D PTV tracking (OpenLPT) and visualization

- OpenLPT: Open-source Lagrangian particle tracking code
 - Based on the Shake-the-Box method
 - Optimized for removing ghost particles at a high particle image density
- Code for trajectory visualization is provided
 - Users can perform further analysis independently



Notes:

- OpenLPT package: [~/Software/3DPTV_Processing_Pipeline/OpenLPT - ShakeTheBox.zip](#)
- OpenLPT post processing: [~/Software/3DPTV_Processing_Pipeline/3D trajectory visualization](#)
- OpenLPT github page: https://github.com/JHU-NI-LAB/OpenLPT_Shake-The-Box
- References: Tan, S., Salibindla, A., Masuk, A. U. M., & Ni, R. (2020). Experiments in Fluids, 61, 1-16.

Date Processing Procedure

Flow Field Imaging Lab




Overview

- Step 1: Extract calibration markers from raw images
- Step 2: Run easyWand to generate the calibration profile
- Step 3: Convert the calibration data to an OpenLPT-compatible format
- Step 4: Use OpenLPT to obtain 3D particle coordinates
- Step 5: 3D trajectory visualization

Step 1: Extract calibration markers from raw images

- Organize raw image files into subfolders by camera under a single calibration folder
 - E.g., “/root_path/calA/cam1”
- Run “droneMarkerTracking_v3.py” in “calibration image pre-processing ” folder
 - GUI pops up and select calibration folder
 - E.g., “/root_path/calA”
 - Images are processed from each camera and the results are displayed
 - White circles and boxes indicate identified markers, while red and green boxes represent the most confident marker identifications. The remaining white circles and boxes are filtered out.
 - Results will be saved in two files
 - wandPts_all.csv (nice to keep)
 - wandPts_clean.csv (cleaned up results with outliers filtered out based on path continuity)



Notes:

- Software path, calibration image pre-processing: ~/Software/3DPTV_Processing_Pipeline/calibration image pre-processing

Step 2: Run easyWand to generate the calibration profile

- Run easyWand5
 - Ensure the software is in “~/Documents/Matlab”.
 - Set this folder as the current directory in Matlab and enter easyWand5 in the command line.
- Compute the initial calibration:
 - Load the previously generated wand points file (wandPts_clean.csv).
 - Click “Load camera profile” and select the appropriate file from the “Camera profiles” folder.
 - Set the wand span to the measured distance between the drone lights.
 - Leave other parameters as default and click “Compute Calibration.”
- Refine calibration
 - Reprojection error will likely still be > 1 pixel, and wand score >>1
 - Press “edit outliers” → pick >one point to be removed with high residual → rerun compute calibration
 - Press “Save results” as .mat file
 - If needed, run Matlab script “cleanEasyWand.m” in “easyWand5/easyWand to OpenLPT” folder
 - load .mat file → filter out points with high error → resave a new calibration .mat file
 - In easyWand5, press “Load previous calibration”, and select newly saved .mat file
 - Press “Compute calibration”
 - Results should yield < 1 pixel reprojection errors → If so, save these results
 - If not, generate a “cleaner” calibration with a lower threshold value



Step 3: Convert the calibration data to an OpenLPT-compatible format

- Create a folder for a new OpenLPT project
 - Add image files, saved in “/cam*/B****.tif” format, where each camera has its own folder
 - Can convert from .raw to .tif using “convertImage.py” in “calibration image pre-processing/Utils”
 - This will be needed to match the config files generated below
- In “easyWand5/easyWand to OpenLPT”, run “easyWand2OpenLPT.m”
 - Will prompt to load the easyWand results file generated
 - This takes the 3D calibration points that easyWand has solved for, and uses them to recompute a new calibration using the Tsai method employed in the OpenLPT calibration tool
 - This also generates the needed configuration files for running OpenLPT
 - Note that scaling factor of 10 is important for OpenLPT not breaking down due to viewing area being too large

Step 4 & 5: Obtain 3D coordinates via OpenLPT and visualize trajectories

- Installation
 - Follow the instructions in the “OpenLPT Instructions from JHU” folder to complete the installation.
- Parameter setup
 - Configuration files are generated automatically, the following two files require particular attention:
 - *trackingconfig.txt*:
 - Keep the view area setting unchanged
 - *iprconfig.txt*:
 - *mindist2d* and *mindist3d*: Set equal to or greater than the calibration error (e.g., ~ 0.25)
 - Set the particle size accordingly.
- Executing the Code
 - Open a command prompt in the OpenLPT root directory: “`~/ShakeTheBox/Release_windows`”
 - Run the executable using the following command with the OpenLPT project root:
 - `.\ShakeTheBox "path-to-OpenLPT-project-root"`
- OpenLPT results
 - Output files are saved as .txt files in the “Tracks” folder
 - Results can be visualized using the scripts provided in the “3D Trajectory Visualization” folder

- The C++ code for OpenLPT needs to be compiled
- In *trackingconfig.txt*: The view area setting should remain unchanged (even though it appears very small, in mm).
- Software path, OpenLPT: `~/Software/OpenLPT – ShakeTheBox`
- Software path, trajectory visualization: `~/Software/3D trajectory visualization`

Demonstration Cases

Flow Field Imaging Lab



Case I: Confetti Settling




Experimental Setup

- This example serves as a concept validation for the 3D PTV system
- Case description and planning
 - Confetti particles made from 6 mm diameter paper hole punches, hand tossed and fall freely
 - Sampling volume: 3 m × 3 m × 2 m, volume center height: 1 m



Raw particle images



Flow Field Imaging Lab




Notes:

- This case was for concept validation in Bristow 2023 EXIF, paper in ~/Field3D_PTV/papers
- Data: ~/Field3D_PTV/Demonstration cases/Case I_Confetti settling
- Reference: Bristow, N., Li, J., Hartford, P., Guala, M., & Hong, J. (2023). Experiments in Fluids, 64(4), 78.


Extract marker coordinates from raw calibration images

- Find “droneMarkerTracking_v3.py” in “calibration image pre-processing” folder
- Select path to raw calibration images and run the code
- Marker coordinates are exported as .csv files for easyWand calibration in the next step

Raw calibration images



Markers identified and coordinates saved as .csv files




Notes:

- Inout: raw marker image; opt put: wand point coordinates
- Color setting: Calibration images are RGB for identifying red and green dots, tracking images are bw
- More details in “Data Processing Procedure”, the code generates .avi file with marker noted for verification
- The calibration working folder, named “calibration\proc\folder,” is included in each demo case.
- Code: ~/Software/3DPTV_Processing_Pipeline/calibration image pre-processing
- Data: ~/Field3D_PTV/Demonstration cases/Case I_Confetti settling/calibration_proc_folder

Compute Calibration Profile

- Run “easyWand5” in Matlab, load wand points and camera profile, enter wand span, compute calibration
- The first result can be with large error, use “Edit outliers” to refine results





Notes:

- Input: wand point coordinates, camera profile, wand span
- Wand span: 34.5 cm
- Camera profile: ~/Software/3DPTV_Processing_Pipeline/Camera profiles
- Code: ~/Software/3DPTV_Processing_Pipeline/easyWand5
- Data: ~/Field3D_PTV/Demonstration cases/Case I_Confetti settling/calibration_proc_folder

Refine Calibration Profile

- Sort calibration points by the residual, remove results that are significantly larger
- Compute calibration again, check wand score and errors, repeat previous step again
- Use “cleanEasyWand.m” in “easyWand5/easyWand to OpenLPT ” folder to filter results if needed
- Save results (.mat file) when the wand score and error meet the desired criteria




Flow Field Imaging Lab



Notes:

- Should target for wand score <1 and pix error <<1
- Output: .mat file from easyWand calibration
- Code: ~/Software/3DPTV_Processing_Pipeline/easyWand5/easyWand to OpenLPT
- Data: ~/Field3D_PTV/Demonstration cases/Case I_Confetti settling/calibration_proc_folder


Calibration Results



- Wand quality score: 0.35
- Wand end point standard deviation: 0.8 mm (total wand length: 34.5 cm)
- Reprojection error of each camera (Cam 1: 0.1 pix; Cam 2: 0.13 pix; Cam 3: 0.12 pix; Cam 4: 0.12 pix)
- Wand score of 1 or less indicates a good calibration (STD/mean = 1%)
- Run code “easyWand2OpenLPT.m” to export config files for OpenLPT

- Input: .mat file from easyWand calibration; Output: config files for OpenLPT
- The code also converts the coordinate system
- Code: ~/Software/3DPTV_Processing_Pipeline/easyWand5/easyWand to OpenLPT
- Data: ~/Field3D_PTV/Demonstration cases/Case I_Confetti settling/calibration_proc_folder

3D Coordinate Calculation using OpenLPT



```
Microsoft Windows [Version 10.0.19045.5854]
(c) Microsoft Corporation. All rights reserved.




C:\Users\Hogan Holography>"C:\Program Files\ShakeTheBox\Release_windows\ShakeTheBox.exe" "E:\Hong_3D_PTV\yue_test\confetti\OpenLPT_proc\trackingconfig.txt"
Code version: W2.1.041619
Reading configuration file...
Select debug mode:
0. No debug
1. Load the 2D position
2. Load the 3D position before shaking
3. Load the 3D position after shaking
4. Load the 3D positions for predictive field
5. Load the predictive field
6. Load the tracks from the initial phase
7. Load the tracks from specific frames
8. Load the tracks from specific frames for back STB
Enter the NO. of the option:0
To save data(1):
E:\Hong_3D_PTV\yue_test\confetti\OpenLPT_proc\Tracks exists
# of outerloops: 4
# of innerloops: 4
IPR on frame 1 of 101
    The 2D particle number in a image is: 969
```

- Create project folder for OpenLPT, organize .tif images into the /cam*/ folders
- Convert easyWand results to OpenLPT configuration files (Use “easyWand2OpenLPT.m” in “easyWand5/easyWand to OpenLPT ”)
- Run OpenLPT (in command window, enter “path-to-ShakeTheBox\ShakeTheBox.exe” “path-to-data\OpenLPT_proc\trackingconfig.txt”)
- Results are saved in “Tracks” folder as txt files

- Make sure path in the configuration files are correct
- The OpenLPT working folder is “OpenLPT_proc_folder”, provided in each demo cases
- Data: ~/Field3D_PTV/Demonstration cases/Case I_Confetti settling/OpenLPT_proc_folder

Trajectory Visualization and Post-processing

- Trajectory yield: 93% of particles successfully captured by the system (464/501)
- Fall speed evaluation: PTV fall speed measurements closely matched independently measured values for individual confetti particles ($\sim 1\text{m/s}$)
- Trajectory analysis: Captured fine-scale motions accurately, reveals details of tumbling motions of the confetti particles



Notes:




- Run `VisualizeResults_2.m` to visualize tracks and analyze velocity distribution
- Location: `~\Code\3DPTV_Processing_Pipeline\3D trajectory visualization\VisualizeResults_3`

Case II: Snow Settling



Experimental Setup

- Investigation on snow settling under different atmospheric conditions
 - Low turbulence (snow morphology play a role)
 - Imaging $4\text{ m} \times 4\text{ m} \times 6\text{ m}$ volume, $\sim 10\text{ m}$ above ground
 - Spatial resolution 6.3 mm/voxel; temporal resolution 200 Hz
 - Duration: in total 26 50-second sequences




Notes:

- This case was the main experiment in Bristow 2023 EXIF, paper in [~/Field3D_PTV/papers](#)
- Data: [~/Field3D_PTV/Demonstration cases/Case II_Snow](#)
- Reference: Bristow, N., Li, J., Hartford, P., Gualala, M., & Hong, J. (2023). Experiments in Fluids, 64(4), 78.

Extract marker coordinates from raw calibration images

- Extract marker coordinates from raw calibration images
 - Find “droneMarkerTracking_v3.py” in “calibration image pre-processing” folder
 - Select path to raw calibration images and run the code
 - Marker coordinates are exported as .csv files for easyWand calibration in the next step




Notes:

- Inout: raw marker image; opt put: wand point coordinates
- Color setting: Calibration images are RGB for identifying red and green dots, tracking images are bw
- More details in “Data Processing Procedure”, the code generates .avi file with marker noted for verification
- The calibration working folder, named “calibration\proc\folder,” is included in each demo case.
- Code: ~/Software/3DPTV_Processing_Pipeline/calibration image pre-processing
- Data: ~/Field3D_PTV/Demonstration cases/Case II_Snow/calibration_proc_folder

Compute Calibration Profile

➤ Compute calibration profile

- Run “easyWand5” in Matlab, load wand points and camera profile (Decimation2X), enter wand spam, compute calibration
- The first result is with large error, use “Edit outliers” to refine results



Notes:



- Input: wand point coordinates, camera profile, wand span
- Wand span: 0.345m
- Camera profile: ~/Software/3DPTV_Processing_Pipeline/Camera profiles
- Code: ~/Software/3DPTV_Processing_Pipeline/easyWand5
- Data: ~/Field3D_PTV/Demonstration cases/Case II_Snow/calibration_proc_folder



Refine Calibration Profile

➤ Refine calibration profile

- Sort calibration points by the residual, remove results that are significantly larger
- Use “cleanEasyWand.m” in “easyWand5/easyWand to OpenLPT ” to filter results
- When desired Wand score is achieved, save results (.mat file) when the wand score and error meet the desired criteria




Notes:

- Should target for wand score ~1 and pix error <1
- Output: .mat file from easyWand calibration
- Code: ~/Software/3DPTV_Processing_Pipeline/easyWand5/easyWand to OpenLPT
- Data: ~/Field3D_PTV/Demonstration cases/Case II_Snow/calibration_proc_folder



Calibration Results

- Calibration result analysis
 - Wand quality score: 1.48
 - Wand end point standard deviation: 3.6 mm (total wand length: 34.5 cm)
 - Reprojection error of each camera
 - Cam 1: 0.33 pix
 - Cam 2: 0.37 pix
 - Cam 3: 0.36 pix
 - Cam 4: 0.34 pix
- Wand score of 1 or less indicates a good calibration (STD/mean = 1%)
- Run code “easyWand2OpenLPT.m” to convert coordinate system and export config files for OpenLPT




- Input: .mat file from easyWand calibration; Output: config files for OpenLPT
- The code also converts the coordinate system
- Every 10th image was included in the data to reduce size. This subset is sufficient for calibration, with minor result variations possible.
- Code: ~/Software/3DPTV_Processing_Pipeline/easyWand5/easyWand to OpenLPT
- Data: ~/Field3D_PTV/Demonstration cases/Case II_Snow settling/calibration_proc_folder

Trajectory Visualization and Post-processing

➤ Process images using OpenLPT

- Create project folder for OpenLPT, organize .tif images into the /cam*/ folders
- Convert easyWand results to OpenLPT configuration files (Use “easyWand2OpenLPT.m” in “easyWand5/easyWand to OpenLPT ”)
- Run OpenLPT (in command window, enter “path-to-ShakeTheBox\ShakeTheBox.exe” “path-to-data\OpenLPT_proc\trackingconfig.txt”)
- Results are saved in “Tracks” folder as .txt files
- Further process the tracking data to extract information; additional details are available in the paper.



Notes:


- Code: ~/Field3D_PTV/Demonstration cases/Case II_Snow settling/3D trajectory visualization
- Data: ~/Field3D_PTV/Demonstration cases/Case II_Snow settling/results_fig
- Reference: Bristow, N., Li, J., Hartford, P., Guala, M., & Hong, J. (2023). Experiments in Fluids, 64(4), 78.

Case III: Pollen Dispersion



Experimental Setup

- Investigation into pollen particle emission from a cottonwood tree
 - Deployments of the system done over Summer 2022 during a time of high pollen density
 - Used 3D PTV system to obtain particle tracks of cottonwood pollen during the daytime
 - A better understanding of pollen emission can potentially improve forecasting for people with allergies



Notes:

- Data: ~/Field3D_PTV/Demonstration cases/Case III_Pollen Dispersion

Extract marker coordinates from raw calibration images

- Extract marker coordinates from raw calibration images
 - Find “droneMarkerTracking_v3.py” in “calibration image pre-processing” folder
 - Select path to raw calibration images and run the code
 - Marker coordinates are exported as .csv files for easyWand calibration in the next step






- Inout: raw marker image; opt put: wand point coordinates
- Color setting: Calibration images are RGB for identifying red and green dots, tracking images are bw
- More details in “Data Processing Procedure”, the code generates .avi file with marker noted for verification
- 1 out of every 10 images was used to reduce the pack size. This reduced set has been shown to be sufficient for identifying the calibration matrix.
- The calibration working folder, named “calibration_proc_folder,” is included in each demo case.
- Code: ~/Software/3DPTV_Processing_Pipeline/calibration image pre-processing
- Data: ~/Field3D_PTV/Demonstration cases/Case III_Pollen Dispersion/calibration_proc_folder

Compute Calibration Profile

➤ Compute calibration profile

- Run “easyWand5” in Matlab, load wand points and camera profile (Decimation2X), enter wand spam, compute calibration
- The first result is with large error, use “Edit outliers” to refine results





Notes:

- Input: wand point coordinates, camera profile, wand span
- Wand span: 0.345m
- Camera profile: ~/Software/3DPTV_Processing_Pipeline/Camera profiles
- Code: ~/Software/3DPTV_Processing_Pipeline/easyWand5
- Data: ~/Field3D_PTV/Demonstration cases/Case III_Pollen Dispersion/calibration_proc_folder


Refine Calibration Profile

➤ Refine calibration profile

- Sort calibration points by the residual, remove results that are significantly larger
- Use “cleanEasyWand.m” in “easyWand5/easyWand to OpenLPT ” to filter results
- When desired Wand score is achieved, save results (.mat file) when the wand score and error meet the desired criteria



Wand score: 0.94459
Wand end point standard deviation: 0.0023042 m
Reprojection errors: 0.20 0.22 0.29 0.20 pixels




Notes:

- Should target for wand score ~1 and pix error <1
- Output: .mat file from easyWand calibration
- Code: ~/Software/3DPTV_Processing_Pipeline/easyWand5/easyWand to OpenLPT
- Data: ~/Field3D_PTV/Demonstration cases/Case III_Pollen Dispersion/calibration_proc_folder



Calibration Results





- Wand quality score: 0.94
- Wand end point standard deviation: 2 mm (total wand length: 34.5 cm)
- Reprojection error of each camera (Cam 1: 0.2 pix; Cam 2: 0.29 pix; Cam 3: 0.29 pix; Cam 4: 0.20 pix)
- Wand score of 1 or less indicates a good calibration (STD/mean = 1%)
- Run code “easyWand2OpenLPT.m” to export config files for OpenLPT

- Input: .mat file from easyWand calibration; Output: config files for OpenLPT
- The code also converts the coordinate system
- Code: ~/Software/3DPTV_Processing_Pipeline/easyWand5/easyWand to OpenLPT
- Data: ~/Field3D_PTV/Demonstration cases/Case III_Pollen Dispersion/calibration_proc_folder

Trajectory Visualization and Post-processing

➤ Process images using OpenLPT

- Create project folder for OpenLPT, organize .tif images into the /cam*/ folders
- Convert easyWand results to OpenLPT configuration files (Use “easyWand2OpenLPT.m” in “easyWand5/easyWand to OpenLPT ”)
- Run OpenLPT (in command window, enter “path-to-ShakeTheBox\ShakeTheBox.exe” “path-to-data\OpenLPT_proc\trackingconfig.txt”)
- Results are saved in “Tracks” folder as .txt files
- Further process the tracking data to extract information; additional details are available in the paper.



Notes:

- Tree shape and location for illustration only
- Due to file size constraints, the demo case includes 200 raw image frames for user testing. Tracking results for these 200 frames are included.
- The full tracking results are also included for visualization demo purposes
- Full tracking results location: `~/Field3D_PTV/Demonstration cases/Case III_Pollen Dispersion/OpenLPT_proc_folder/Tracks_Full`
- Data: `~/Field3D_PTV/Demonstration cases/Case III_Pollen Dispersion/3D trajectory visualization`